



KangooFix™

Neonatal Restraint System



The **KangooFix Neonatal Restraint System** was developed to safely and effectively cradle and secure a newborn during non-critical ambulance transportation, allowing the baby to travel with their parent.

Not only does this eliminate the need to use two vehicles and crews for transfers, but it also ensures a safer and improved environment for both parent and baby throughout the journey.

- **EN 1789 compliant** for safety
- **Skull cap included** to help maintain baby's body heat
- **Quick release system** allows fast access for EMS professionals
- **"Face-up" mode** enables EMS professionals to quickly monitor newborn



Safely secures newborn to parent, maintaining close contact at all times



Inner harness and warming cover ensure breathability and comfort



Material was clinically tested to guarantee no negative reaction on sensitive baby skin

KangooFix

Specifications



ITEM NUMBER 0374906

LENGTH (MAX) 43.3 in | 110 cm

LENGTH (MIN) 9.8 in | 25 cm

WIDTH (MAX) 36.6 in | 93 cm

WIDTH (MIN) 5 in | 13 cm

HEIGHT (MAX) 1 in | 2.4 cm

HEIGHT (MIN) 1.6 in | .4 cm

WEIGHT 2.2 lb | 1 kg

WEIGHT CAPACITY (MAX) 11.1 lb | 5 kg

WEIGHT CAPACITY (MIN) 3.5 lb | 1.6 kg

MATERIAL "Breath-o-prene" Oeko-Tex 100, polyester, metal, and plastic

Features

- **5-point outer harness** attaches to the ambulance cot harness
- **Inner harness and heat cover** are sewn in to enable secure fit
- **Snug cocoon design** and integrated skull cap keeps baby warm and secure
- **Close contact** between parent and baby is clinically proven to improve condition of both parent and baby during transfers
- **Lightweight** and easy to store
- **Easy to clean**

Compliance

The **KangooFix** has been designed to meet international test standards to increase safety for you and your patient.

- EN 1789