



BABYPOD 20 Infant Transport System™



- **Greater visibility and access baby** – improved larger viewing window area and easier to access
- **Smarter sliding top-lid** – initially slides back halfway to allow easy access, then opens fully by lifting the lid upwards, designed for use in confined & narrow spaces. Equipment such as ventilators also can be positioned at the side of the Pod without hindering access
- **CO2 evacuation improves baby's air supply quality** – an external inlet point on the head-end of the Pod allows air or oxygen mix to be fed into the unit, thereby diminishing any build-up of CO2
- **Better access for clinical intervention** – larger access areas for tubing & cables now allow more medical equipment to be used simultaneously
- **Revised trolley strapping system** – uses high-tensile webbing and metal quick-connect buckles all tested to 9kN, suitable for both air and road transport. Straps are now attached to and accessible from the outside of the Pod, making cleaning quicker and easier
- **Safer snug cocoon area for precious cargo** – vacuum mattress has slots in edges to allow the internal positioning straps to be fed through the mattress, thereby allowing the straps to better secure the mattress around the baby.

BabyPod 20 incorporates new technological developments and additional user-requested features. Bringing to market an extra safe high-tech infant transport device to meet the time critical challenges and treatment requirements usually associated with poorly babies.

Following on from it's successful predecessor BabyPod II EV04, the redesigned BabyPod20's most striking difference is the completely re-engineered slide-lift lid, giving easier access in narrow confined spaces.

The viewing area in the BabyPod20 is larger and easier to access than ever before, giving better visibility and access to the patient, particularly in restricted environments such as the interior of rotary or fixed wing aircraft where space is at a premium.



BABYPOD 20

Infant Transport System™



Features & Benefits

- **Simple to transport** – easily lifted in and out of vehicle, optional carry bag also available
- **Unrivalled stability** – using vacuum mattress and inner patient-positioning safety straps
- **Exceptional versatility** – ideal for use in hospital transfers, relief operations, emergency midwifery call-outs, neonatal revival cases, in ambulance, car or air transfers
- **Convenient to use** – needs no electricity, maintain temperature by use of TransWarmer heat pads
- **Proven performance** – crash tested to EN 1789, CEN compliant. Fixing restraints tested to 9kN for use in aircraft
- **Highly compatible** – will fit on all models of ambulance trolley, can be used with all major ventilators enabling the addition of improved O2 environment where required
- **Cost-Saving** – at a cost that is less than 20% of a standard transport incubator, with reduced oxygen consumption due to its compact interior and lighter due to its carbon fibre construction

Specifications

LENGTH	1000 mm
WIDTH	480 mm
HEIGHT	320 mm
WEIGHT	8,6 kg
incl. vacuum mattress	9,6 kg
CAPACITY	2 - 8 kg / 4 - 17 lbs

Made of Carbon fibre
Four Layers of Protection
Carbon Fibre (F1) Shell
Shock-absorbing laminated foam interior
Unique design Vacuum positioning mattress
Vibration absorbing Gel warming mattress

Approvals and compliance:
CE and FDA approved
CEN compliant, dynamically crash tested to BS EN 1789
ISO 9001:2000 Quality Management
ISO 13485:2016 Design + Manufacturing of Medical Devices

Accessories/Spares



- Replacement Vacuum Mattress BP0D2-002



- Additional positioning Straps (single use) BP0D2-003



- Replacement Shock Absorbent Foam BP20-005



- Additional 9kN Trolley Fixing Straps BP20-006



- Replacement Lid Assembly BP20-007



- Protective Carry / Storage Bag BP20-013



- Vacuum Pump for Vacuum Mattress BP0D2-020

Part Numbers

BP20-001	BabyPod 20 (includes Vacuum Mattress, Shock Absorbent Foam Set, 1 set Positioning Straps (single use), Stretcher Fixing Straps, Protective Dust Cover & 1 TransWarmer)
BP20-013	Protective Storage/Carry Bag for BabyPod 20
BP0D2-020	Hand Vacuum Pump for vacuum mattress (for BabyPod 2 & BabyPod 20)
BP0D2-002	Replacement Vacuum Mattress (for BabyPod 2 & BabyPod 20)
BP0D2-003	Additional single-use Infant Positioning Strap Set (for BabyPod 2 & BabyPod 20)
BP20-005	Replacement Shock Absorbent Inner Foam Set for BabyPod 20 only
BP20-006	Trolley Fixing Strap Set for BabyPod 20, metal buckles, tested to 9kN
BP20-012	Replacement Protective Dust Cover for BabyPod 20
BP20-007	Replacement Lid for BabyPod 20
TRANS06	Additional TransWarmer heating pads (pack of 6)