

FERNO EasyFix PLUS II Vacuum Mattress

Upgraded design for improved product durability, providing further enhanced immobilisation, eliminating the need to use a cervical collar, and minimised storage space



**FOR MORE DETAILED
INFORMATION SEE OVER**

Specifications

Length	Width (max)	Thickness	Weight	Load Limit
2080 mm	720 mm 1300 mm	45mm	8.2 kg (incl straps)	250 kg/39st/550lbs

FERNO – is a world leading innovator and manufacturer of patient transportation solutions.

Our product range is born of a closer understanding of our customers' needs and the practical demands of their work. Ferno's industry-leading **EasyFix PLUS II Vacuum Mattress** typifies our passion for top quality, highly functional design... patient transfer products that don't just do the job, they do it **BETTER**.

call +44 (0) 1274 851 999

email sales@ferno.co.uk

www ferno.co.uk



FERNO EasyFix PLUS II Vacuum Mattress

Following extensive user feedback we have further upgraded and enhanced the EasyFix PLUS to deliver even better head, neck and upper back immobilisation and enhanced durability of product when using in harsh operating situations.

FERNO's EasyFix PLUS II provides an effective and secure vacuumed environment for a safer transfer of patients and casualties.

Designed to give **better product protection in unforgiving operating environments**, the ground-facing base surface made from a durable low-friction glass-fibre material now runs along the full length of the vacuum mattress, making it even easier to manoeuvre and gives even more protection when dragging or bringing someone downstairs.

The V-shaped inner chamber system and rigid base gives effective **C-spine immobilisation without the need to use a cervical collar**, which can cause extra trauma to the patient during positioning of a collar (can also be used with a collar).

Improved design of the inner head/neck support board gives better protection when moving on the ground or down stairs and reduces wear over time.

An **enhanced angled headstrap** on the EasyFix II models also makes deployment even easier and more effective.

The **EasyFix PLUS II** vacuum mattress raises the level in patient immobilisation & transportation equipment, allowing a patient to be moved from scene all the way through to scanner without the need to remove then when transferring between transportation devices, in all environments and situations.

- **Improved design** in the head and neck area for superior immobilisation of patient, with no collar required
- **Durable and highly effective** – hard-wearing base fibreglass material reduces friction, making surface movement of the vacuum mattress easier
- **Integrated immobilisation** – integral head straps, X-Strap and across-the-shoulder flaps to hold head in place
- **High carrying capacity** – load limit 250kg
- **X-ray transparent, MRI compatible** – no need to move the patient out of the vacuum mattress for scanning
- **Easy storage** – folds easily to take up 20% less space than other designs of vacuum mattress

*Standard Green/Yellow FERNO EasyFix II also available



Part Nos:

271511101	EasyFix PLUS II Vacuum Mattress "Ready2Go kit" (incl. pump, repair kit & carry bag)
271512101	EasyFix PLUS II Vacuum Mattress "BackPack Ready2Go kit" (incl. pump, repair kit & backpack)
271511100	EasyFix PLUS II Vacuum Mattress (incl. restraints only, no pump or bag)
271000708	Backpack for EasyFix and EasyFIX PLUS range (Mkl & II), Red



FERNO - EasyFix PLUS II Vacuum Mattress FPC001 03.21

call +44 (0) 1274 851 999
email sales@ferno.co.uk
www ferno.co.uk

