

FERNO EasyFix Vacuum Mattress

Revolutionary new body shape sets new standards in vacuum mattress immobilisation



**FOR MORE DETAILED
INFORMATION SEE OVER**

Specifications

Length	Width	Weight	Load Limit
2080 mm	1300 mm (Max)	8.5 kg (incl straps) 9.6 kg (incl straps, footrest & pump)	250 kg/39st/550lbs

FERNO – is a world leading innovator and manufacturer of patient transportation solutions.

Our product range is born of a closer understanding of our customers' needs and the practical demands of their work. Ferno's industry-leading **EasyFix Vacuum Mattress** typifies our passion for top quality, highly functional design... patient transfer products that don't just do the job, they do it **BETTER**.

call +44 (0) 1274 851 999
email sales@ferno.co.uk
www ferno.co.uk

 **FERNO**[®]

FERNO EasyFix Vacuum Mattress

Achieve secure immobilisation with improved vacuum mattress design

EasyFix utilises advanced technological developments to deliver superior immobilisation, exactly where and when it's needed. The system uses a new V-shaped body design to optimise the vacuum mattress capabilities, whilst at the same time drawing on enhanced internal chamber technology for superior stability throughout. The revolutionary X-shape restraint system immobilises the upper body and/or hip and pelvis, demonstrating the versatility of this unique mattress.

However sophisticated an immobilisation technique may be, it also needs to be quick and easy to apply. The **EasyFix** offers improved performance here too, as the suction devices fits directly onto the mattress valve so no time is wasted attaching adaptors – and easily adjusted strap mechanisms ensure operators can save time when it's needed most.

- **Superior immobilisation** – unique chamber design and body shape with built-in wooden upper supports to hold the patient in the correct position more securely, head and shoulder support also provides further patient safety
- **Improved immobilisation of upper back, neck and head** – reinforced base plate
- **Durable and highly effective** – double layered material helps maintain vacuum longer
- **High carrying capacity** – load limit 250kg
- **X-ray and MRI transparent** – no need to move the patient to another device for diagnostics
- **Easy storage** – folds easily to take up 20% less space than other designs of vacuum mattress
- **Simple maintenance** – polyester fibre, PVC coated
- **Smooth handling** – with 10 ergonomically designed handles
- **Enhanced use** – X-Ray and MRI translucent



Part Nos:

271501000 EasyFix Vacuum Mattress

271501001 EasyFix Vacuum Mattress Kit (includes footrest pump, repair and bag)

FERNO - EasyFix Vacuum Mattress FPC003 10.15

call +44 (0) 1274 851 999
email sales@ferno.co.uk
www ferno.co.uk



Ferno (UK) Limited, Stubs Beck Lane, Cleckheaton, West Yorkshire, BD19 4TZ England.