

# **FERNO PXB**

## **Pool Extraction Board**

The faster and safer extrication standard for swimming pool Lifeguarding that requires only 2 operators



Developed by **FERNO**  
in partnership with

**ROYAL  
LIFE SAVING  
SOCIETY UK**  
[www.rlss.org.uk](http://www.rlss.org.uk)

**FOR MORE DETAILED  
INFORMATION SEE OVER**

### **Specifications**

Length	Width	Depth	Weight	Load Limit
2000 mm	400 mm	40 mm 200 mm with head blocks	11 kg	182 kg (28stone)

**FERNO** – is a world leading innovator and manufacturer of patient transportation solutions.

Our product range is born of a closer understanding of our customers' needs and the practical demands of their work. Ferno's industry-leading **PXB Pool Extraction Board** typifies our passion for top quality, highly functional design... patient transfer products that don't just do the job, they do it **BETTER**.

**call** +44 (0) 1274 851999  
**email** [sales.uk@ferno.com](mailto:sales.uk@ferno.com)  
**www** [fernouk.com](http://fernouk.com)



# FERNO PXB Pool Extraction Board

Now lifeguards can perform a significantly faster but controlled shallow or deep-water, deck-level or free-board pool rescue, with reduced injury risk to themselves and improved safety for the casualty – with two-person operation in a deck-level pool and recommended three-person operation in a free-board pool (not exceeding a height of 30 cm).

## The PXB is packed full of innovations

The neutral-buoyancy of the board is such that a casualty can be held on the water surface without the need for any floatation devices. The restraint straps are made from durable polypropylene material, incorporating an easy-to-operate fixation system using specially developed anti-bacterial VELCRO®, the headblocks & detachable headblock baseplate also containing anti-bacterial additives for superlative infection control.

Durable edge-stitched Neoprene anti-slip mat is specifically designed for use in water and the rubberised anchor mat ensures safe and effective stabilisation of the board to the pool-side when in operation. The board underside is designed to give a firm anchor against poolside whilst casualty is being brought onto the board and restrained/immobilised, and then dragged out of the pool side using the attached anchor mat.

## With many additional benefits

- Fewer lifeguards required
- Easier to rescue casualty from pool
- Provides full C-spine Immobilisation
- Faster compared to other systems
- Simpler and more intuitive to use
- Restraints colour-coded for ease of deployment
- Reduces manual handling risks
- Less training required
- Easier to maintain and keep in an operational condition



PXB detachable headpad with Velcro pad for holding headblocks in place

### Part Numbers

010776700 Ferno PXB board with detachable headblock baseplate, including wall bracket

### Accessories and replacement parts

038245650	Replacement headblock baseplate
038245300	Replacement set of 2 headblocks & 1 headstrap
038245250	Replacement head-strap only
038244500	Full set of PXB restraint straps (5 in total – green, blue, yellow, black & red)
038245350	Replacement PXB Chest/Waist strap, Green
038245400	Replacement PXB Chest/Waist strap, Blue
038245450	Replacement PXB Leg strap, Yellow
038245500	Replacement PXB Leg strap, Black
038245550	Replacement PXB Ankle strap, Red
0761-3001	Replacement PXB neoprene anti-slip mat
0761-0119	PXB board wall-mounted fixing bracket
0761-0102	Replacement anchor mat

ALL FERNO products come with a comprehensive 'User Manual', for guidance on the recommended use of the product, However accredited training must be undertaken. Ferno (UK) Ltd cannot be held responsible for equipment that has not been serviced and maintained in accordance with its guidelines, specification or for the incorrect use of equipment under working conditions.

FERNO - PXB MkII FPC001 08.25

**call** +44 (0) 1274 851999  
**email** sales.uk@ferno.com  
**www** ferno.com



Ferno (UK) Limited, Stubbs Beck Lane, Cleckheaton, West Yorkshire, BD19 4TZ, England