

FERNO Compact 2 Track



In-vehicle storage bracket

Cabin-mounted two-part bracket holds both the chair and track - so they're to hand when you need them and secure when you don't



FOR MORE DETAILED
INFORMATION SEE OVER

Specifications

FERNO – is a world leading innovator and manufacturer of patient transport solutions.

Our product range is born of a closer understanding of our customers' needs and the practical demands of their work. Ferno's industry-leading **Compact 2 Track Bracket** typifies our passion for top quality, highly functional design... patient transfer products that don't just do the job, they do it **BETTER**.

Width	Depth with Chair	Height Empty	Height with Chair	Total Height Required	Mass without Chair
435 mm	250 mm	673 mm	800 mm	850 mm for safe removal	7.9kg

call +44 (0) 1274 851 999
email sales@ferno.co.uk
www ferno.co.uk



FERNO Compact 2 Track In-vehicle storage bracket




Secured, ready for use and easy to deploy

If you carry the Compact 2 Track in your Ambulance this crash-tested two-piece storage bracket is an essential extra.

The track fits in first and can stay in the bracket if you need to use the Compact 2 purely as a carry chair, the chair locks in over the track and is secured further with a strap. CEN Crash tested to 10G and CEN compliant, you can rest assured that when stored, the equipment is safe and secured to the highest standard.

- **Purpose** – designed to secure the Compact 2 Track in vehicle
- **Protects** – tested and certified to BS EN 1789
- **Easy to use** – the track and chair simply click into place and are easily deployed
- **Easy to clean and maintain** – ask about our comprehensive service contract for all equipment



Part Nos:

009103800

Compact 2 Track 2-Piece In-Vehicle Storage Bracket

FERNO - C2T Bracket FPC008 11.17

call +44 (0) 1274 851 999
email sales@ferno.co.uk
www ferno.co.uk



Ferno (UK) Limited, Ferno House, Stubs Beck Lane, Cleckheaton, West Yorkshire, BD19 4TZ England.