

FERNO MEGASUS Trolley

The trolley that eliminates the need for a specialised vehicle to attend bariatric patients



FOR MORE DETAILED
INFORMATION SEE OVER

Specifications

Length	Width	Height	Weight	Load Limit
1950 mm 1680 mm (With folded head)	570 mm 955 mm (Extended)	880 mm (Raised) 460 mm (Lowered)	73 kg	200kg/31st/434lbs 400kg/63st/882lbs (Lowest Position)

FERNO – is a world leading innovator and manufacturer of patient transportation solutions.

Our product range is born of a closer understanding of our customers' needs and the practical demands of their work. **The Megasus Trolley** typifies our passion for top quality, highly functional design... patient transfer products that don't just do the job, they do it **BETTER**.

call +44 (0) 1274 851 999
email sales@ferno.co.uk
www ferno.co.uk

 **FERNO**

FERNO MEGASUS Trolley

Transport larger patients safely with better comfort to the patient



Bariatric trolleys need to perform two vital functions – keep the patient comfortable and secure, and ensure that the operator can perform their job easily and without risk of injury. **The Megasus** is Ferno's highly effective solution to this challenge, designed to carry patients up to 400kg and accommodate their associated body mass.

The trolley's configurability, high load capacity and excellent manoeuvrability also afford operators another major advantage – it allows a standard ambulance to attend call-outs which might otherwise have required a specialised vehicle.

- **Lift heavier patients smoothly and easily** – 400 kg in lowest position, 200 kg raised
- **Additional comfort and space** – extendable cotsides, allowing a larger mass capacity, as and when required
- **Excellent security** – Complete four-point body harness and leg restraints
- **Operation without strain** – dual action raise/lower hydraulics, gas-strut assisted backrest, extendable push/pull handles at both ends, optional detachable push-poles at foot-end



Megasus with cotsides fully extended



Megasus with option of detachable push-poles at foot-end

Part Nos:

MEGBIO-STD-AB/PR	Megasus c/w ABPR Mattress, extending cotsides, Biosafe restraints, Standard 2 Part Lock
MEGBIO-LHS-AB/PR	Megasus c/w ABPR Mattress, extending cotsides, Biosafe restraints, LHS-operated 2 Part Lock
MEGBIO-NMI-AB/PR	Megasus c/w ABPR Mattress, extending cotsides, Biosafe restraints, 2 Part NMI Track Lock
MEGBIO-SNMI-AB/PR	Megasus c/w ABPR Mattress, extending cotsides, Biosafe restraints, LHS 2 Part NMI Track Lock
MEGBIO-UNW-AB/PR	Megasus c/w ABPR Mattress, extending cotsides, Biosafe restraints, 2 Part Unwin Track Lock
MEGBIO-LHUNW-AB/PR	Megasus c/w ABPR Mattress, extending cotsides, Biosafe restraints, LHS Unwin Track Lock
MEGBIO-SL150-AB/PR	Megasus c/w ABPR Mattress, extending cotsides, Biosafe restraints, SlideLock (150mm)
MEGBIO-SL255-AB/PR	Megasus c/w ABPR Mattress, extending cotsides, Biosafe restraints, SlideLock (255mm)

Also available with webbing restraints instead of Biosafe, p/n starting 'MEG-' in place of 'MEGBIO-'

Folding push-poles at foot end also available as an in-built option. Request for details.

SlideLock also available with bespoke distance between lock positions according to vehicle size/requirements.

Options & Accessories:

009001700	Detachable IV Pole for Megasus
008551000	Oxyclip cylinder holder for CD cylinder
009003410	Base mounted dual-strap cylinder holder
082-2097	Patient shield for trolley
031346990	Trolley dust cover, grey AB material

FERNO - MEGASUS FPC006 02.21

call +44 (0) 1274 851 999

email sales@ferno.co.uk

www ferno.co.uk



Ferno (UK) Limited, Ferno House, Stubs Beck Lane, Cleckheaton, West Yorkshire, BD19 4TZ England.