

NAJO[®]

BACKBOARDS



Najo Backboards provide medics and rescue personnel with the versatility to respond to any situation.

Made from a unique blend of advanced plastics and foams, Najo boards are strong, yet lightweight and meet OSHA recommendations. The one-piece, HDPE construction is impervious to all bodily fluids, easy to decontaminate, and 100% radiolucent.

Najo boards are buoyant and can float a patient in water



The **Najo[®] Lite Board** is the lightest Najo backboard at just 14.5 lbs.



The **Najo[®] RediHold** features two handles at tapered end to allow for more stable transport

- **Guaranteed for life under normal usage**
(except for the Najo Lite Board, which comes with a two-year warranty)
- **Angled edges** allow for smooth log-rolls
- **100% Radiolucent**
- **Meets OSHA recommendations**

Specifications



Orange No Pin
No. 275201108

Orange 10 Pin
No. 275201208



Yellow No Pin
No. 275201105

Yellow 10 Pin
No. 275201205



Orange No Pin
No. 275101108

Orange 10 Pin
No. 275101208



Yellow No Pin
No. 275101105

Yellow 10 Pin
No. 275101205

	Najo Lite Board™	Najo RediHold®
LENGTH (MAX)	72 in 1830 mm	72 in 1830 mm
WIDTH	16 in 406 mm	16 in 406 mm
THICKNESS	1.75 in 44 mm	1.75 in 44 mm
WEIGHT	14.5 lb 6.6 kg	16 lb 7.2 kg
LOAD LIMIT	450 lb 204 kg	600 lb 272 kg
FLOATATION AT 1.2 BUOYANCY	154 lb 69 kg	136 lb 62 kg

Features

Najo Lite Board

- Lightest Najo board at just 14.5 lb
- Central slot in low half for immobilizing each leg separately
- 2-year warranty

Najo RediHold

- Two handles at tapered end for more stable transport
- Smooth head area to insure that all head immobilization devices adhere
- Guaranteed for life under normal usage

Accessories



Universal Head Immobilizer
No. 0313855



Wizloc™ Cervical Collar
No. 0822074 (3 pack)



HeadHugger™ Disposable Head Immobilizer
No. 0832092 (single)
No. 0832093 (10 pack)



Fastrap™ Restraint System
No. 0313916



Model 430 Restraints
Please contact your local sales representative for product details



Model 430P Restraints
Please contact your local sales representative for product details



Model 439 Restraints
Please contact your local sales representative for product details



Model 436 Speed Clip Restraints
No. 0314118 (set of 4, 5' w/metal buckles)
No. 0314119 (set of 4, 7' w/metal buckles)
No. 0314120 (set of 4, 5' w/plastic buckles)
No. 0314121 (set of 4, 7' w/plastic buckles)