

# **FERNO Compact 5 Track**

Designed with both patient safety and comfort in mind, this multi-purpose chair is suitable for both evacuation and regular daily PTS use



**FOR MORE DETAILED  
INFORMATION SEE OVER**

## Specifications

**FERNO** – is a world leading innovator and manufacturer of patient transport solutions.

Our product range is born of a closer understanding of our customers' needs and the practical demands of their work.

Ferno's industry-leading **Compact 5 Track Chair** typifies our passion for top quality, highly functional design... patient transfer products that don't just do the job, they do it **BETTER**.

Height	Depth	Width	Weight	Load Limit
1160 mm (folded) 510 mm (seat height)	300 mm (folded) 930 mm (track deployed)	515 mm (folded)	17kg (chair+track)	200kg/31st/434lbs
Track module length	Track module width	Track module depth		
730 mm	340 mm	120 mm		

**call** +44 (0) 1274 851 999  
**email** sales@ferno.co.uk  
**www** ferno.co.uk



# FERNO Compact 5 Track

## Multi-function track and carry chair, designed for a wide range of applications and uses

The Compact 5 Track is based on the same proven design of its predecessor, but improved to give better patient safety and comfort, whilst also providing enhanced operator functionality. The built-in tracks allows a patient to be moved down stairs with ease, the chair's 4 wheels making horizontal movement effortless.

Robust design for high usage makes this the perfect chair for Patient Transport Services patient transfers, as well as for use as an evacuation chair.

- **Multi-functional** - designed to be used as a carry chair, tracked stair chair and wheelchair
- **Enhanced Patient Comfort** - includes hinged armrests and folding footrest
- **Added Security** - chest, waist, and ankle biosafe patient restraints
- **Handles brilliantly in compact spaces** - front castor wheels and rear braked wheels for ease of manoeuvrability
- **Outstanding versatility** - retractable lower grab handles and adjustable rear lifting handle to accommodate different height operators



### Part Number

071481701	Compact 5 Track Chair
081931500	Wall Bracket
038299801	Dust Cover

FERNO - Compact 5 Track FPC005 02.16

**call** +44 (0) 1274 851 999  
**email** sales@ferno.co.uk  
**www** ferno.co.uk



Ferno (UK) Limited, Stubbs Beck Lane, Cleckheaton, West Yorkshire, BD19 4TZ England.