

# COMPACT POWERTRAXX

*Affordable Powered Track Chair for everyday use*



## **Compact PowerTraxx Chair - bringing POWER to day-to day patient transfers**

One of the most popular Ambulance Chairs used in the UK is the FERNO Compact 2 Track manual track chair, widely used every day in many situations. It became the go to solution to get a patient downstairs and into the back of the Ambulance. Lightweight and compact in size, it is extremely manoeuvrable around tight bends and when navigating narrow stairways, reliably used both in Frontline emergency service patient transfers and PTS planned transfers and patient repatriations.

But we all know that lifting when moving a patient **UP or DOWN** stairs has its risks. **NOW YOU HAVE THE POWER** to eliminate lifting with the **Ferno Compact PowerTraxx Powered Track chair**.



**Compact PowerTraxx** - a robust, lightweight and compact Powered Track Chair for transferring patients both up and down stairs, eliminating the need to lift, giving power-assistance to crews when moving up ramps and inclines, also contributing to an overall reduction in AMBULANCE vehicle weight.

## Features

- Can use as a standard carry chair as well as a Powered Track chair
- Attachable Powered Track module
- Control box connects wirelessly to motor module
- Variable and fixed-speed controls
- Feature allows chair to be moved downstairs even if run out of power
- Light to illuminate dimly-lit stairways
- Adjustable-height rear handle
- Front operator grip handles
- Removable wipe-clean seat cover
- Patient head pad
- Large size wheels
- Biosafe chest & ankle patient safety straps

## With many Benefits

- Not only **POWER UP and DOWN** stairways, you can now **LOAD and UNLOAD** the patient into a vehicle using the Ambulance Ramp or rear/side steps under **POWER**, and also use the **POWER** to assist on inclines, curbs, doorsteps, and other such obstacles.
- The **COMPACT PowerTraxx** chair can be used on its own as a standard carry chair for "non-lift" wheeled transfers, simply attaching the PowerTraxx module when required and eliminating the need to lift
- Only 23.8 kg, working hard to reduce the overall equipment weight within the AMBULANCE
- Easy use touch **CONTROL** technology gives you **PROVEN LITHIUM-ION BATTERY POWER**, how you want it and when
- Training is kept to a minimum where crews have already been trained and fully familiar with the Compact 2 Track chair
- Same controls and similar in operation to the Venice PowerTraxx makes transition to Compact PowerTraxx easy
- Proven design, being exceptionally manoeuvrable in tight spaces and difficult stairways, has made the Compact range the most used Ambulance chairs in the UK and now is available in a **POWERED** version.....



## Specifications

### COMPACT PowerTraxx

CHAIR HEIGHT FOLDED (without track attached)	850 mm
CHAIR DEPTH FOLDED (without track)	290 mm
CHAIR HEIGHT UNFOLDED	1102 mm
CHAIR HEIGHT WITH REAR HANDLE FULLY EXTENDED	1485 mm
SEAT HEIGHT	485 mm
CHAIR WIDTH	450 mm
TRACK ASSEMBLY LENGTH	720 mm
TRACK ASSEMBLY WIDTH	340 mm
WEIGHT OF CHAIR FRAME	13.3 kg
WEIGHT OF POWERED TRACK MODULE (INCL BATTERY)	10.0 kg
TOTAL WEIGHT	23.3 kg
SWL / Max load capacity (31.5 stone)	200 kg
Ferno Compact PowerTraxx in-vehicle storage bracket	3.9 kg

### Standards:

- Compliance to EN 1865:2015
- Compliance to Medical Device Regulation (EU) MDR 17/745

## Part Numbers

**071481801** Compact PowerTraxx chair c/w 1 battery

## Accessories



Spare 5A Milwaukee M28 Battery

75-0710-001



2 part In-Vehicle Storage Bracket

009103900