

# FERNO

## S-CAPEPLUS® Inflatable

Evacuation device for where storage space is limited

For the fast but safe movement of casualties and patients in difficult situations



FOR MORE DETAILED  
INFORMATION SEE OVER

### Specifications

**FERNO** – is a world leading innovator and manufacturer of patient transport solutions

Our product range is born of a closer understanding of our customers' needs and the practical demands of their work.

Ferno's industry-leading **S-CAPEPOD** and **S-CAPEPLUS** evacuation devices typifies our passion for top quality, highly functional design patient transfer products that don't just do the job, they do it **BETTER**.

S-CAPEPLUS	In carry bag	Base size	Dims Flaps Open	Weight	Load Limit
S-CAPEPLUS Standard	600x350x600 mm	2100*600 mm	1900x2850 mm	5.5 kg	200 kg
S-CAPEPLUS Inflatable	TBA	TBA	TBA	TBA	TBA

**call** +44 (0) 1274 851 999  
**email** sales@ferno.co.uk  
**www** ferno.co.uk



# FERNO

## S-CAPEPLUS® Inflatable

For the fast but safe movement of patients and casualties alike from difficult locations and situations.

Based on the same concept as S-CAPEPLUS Standard, **S-CAPEPLUS Inflatable** is designed to minimise the space required to store the evacuation device. With an inflatable base instead of a foam mattress, the S-CAPEPLUS Inflatable is compact in size, easy to carry and takes up little space when stored.

Principally designed for use in the Ambulance sector, but also applicable for use by all Emergency Services organisations for moving people out from difficult locations. S-CAPEPLUS Inflatable can also be used as a general evacuation device where a chair or stretcher cannot be used, complimentary to the range of evacuation equipment available to use, thereby covering all eventualities and situations.

The unique design of using flaps of material at each side and at the foot-end of the inflatable base, which fold over each other to create a sturdy 'bag' in which to envelop the occupant. Using a heavy-duty hook & loop fixation system to hold it all together, it is quick and easy to deploy, giving an extremely fast, safe and effective system of evacuation or movement of a patient or casualty in all scenarios.

Hand-hold straps at the head and foot-end allow the occupant to be easily and safely manoeuvred along the floor and down stairs. Made with a special high-strength material base, the S-CAPEPLUS Inflatable glides easily during evacuation, but giving sufficient grip to ensure a controlled descent when moving down a set of stairs or slope.

Combined with Ferno's range of patient carry chairs, manual and powered track chairs, it ensures that people of all sizes and conditions can be effectively moved in any emergency situation, seated or recumbent.

- **Compact in size** - inflatable base taking up little space when stored
- **Fast and safe** – easy to deploy, quick to move patients under control with no lifting
- **Extremely easy to use** – two-person operation is advised but in extreme emergency can be operated by one person only
- **Compact in size** - inflatable base taking up less space when stored
- **Versatile** - can be used in a wide range of scenarios
- **Long service life, low maintenance** – strong and hard-wearing, easy to clean
- **Infection control** - made from materials with AB additives, cover removable for high-temperature washing



Compact in size, easy to store and carry



Folding side flaps for securely enveloping the patient



Inflatable base



Foot section for support and comfort



Head-end & foot-end H-straps for controlled movement of patient



Internal restraints for securing the patient

FERNO - S-CAPE INFLATABLE FPC001 04.20

#### Part Number

5NINFL6021 'S-CAPEPLUS Inflatable' evacuation device with carry/storage bag and pump

call +44 (0) 1274 851 999 email [sales@ferno.co.uk](mailto:sales@ferno.co.uk)  
www [ferno.co.uk](http://ferno.co.uk)

Ferno (UK) Limited, Stubs Beck Lane, Cleckheaton, West Yorkshire, BD19 4TZ England.

