

# **TRANSCEND™**

EVACUATION CHAIR

## A more comfortable ride.

The Transcend™ Evacuation Chair helps you quickly and safely transport a seated passenger down stairs and over flat surfaces in an evacuation situation.



PROVIDES COMFORTABLE RIDE ON ALMOST ANY SURFACE WITH 60% LESS SHOCK THAN COMPETITIVE STAIR CHAIR

### EXTENDING FOOTREST

OFFERING ADDITIONAL SUPPORT AND COMFORT FOR PATIENTS WHEN EITHER STATIONARY OR MOVING.

### REAR LIFT HANDLES

WITH IMPROVED ERGONOMICS FOR LIFTING OR CARRYING THE EVACUATION CHAIR



## ENHANCED PATIENT CARE, COMFORT, AND SAFETY

- Extending footrest and lower track angle increase patient sense of security
- Patient seating surface is contoured and has removable injection molded seat panels
- Patient restraints for security and safety
- Armrests offer additional support and security for patients while enhancing comfort and positioning during transport or transfers\*
- Headrest provides added comfort and stability for the patient's head and neck during transport, helping reduce movement and increase confidence\*

## INCREASED PARAMEDIC SAFETY, EFFICIENCY, AND MANEUVERABILITY

- Smooth 5-inch front swivel wheels and extra tall 8-inch rear wheels with integrated shock absorbing tires that navigate rough terrain and do not collect debris
- 5-position lift bar and handles reduce the strain of lifting and provide maximum control for operators of various heights (wireless motor controls work in any handle position)\*
- One-inch track-to-ground clearance provides smooth rolling over carpet and rough surfaces, as well as easier descent down stairs
- Rear lift handles provide improved ergonomics for safer, more comfortable lifting and greater control\*
- Positive locking hinge mechanism ensures chair is locked in place whether deployed for use or folded for carrying and storage

## ENGINEERED FOR RELIABILITY, DURABILITY, AND PERFORMANCE

- All-metal frame construction provides durability and long life
- 15 year warranty

\*OPTIONAL ACCESSORIES

## SPECIFICATIONS

ITEM NUMBER\*  
**0731404**

HEIGHT (MAX)  
**64" / 1633 MM**

HEIGHT (MIN)  
**40" / 1022 MM**

WIDTH (OVERALL)  
**20" / 505 MM**

WIDTH (SEAT)  
**16" / 408 MM**

DEPTH (FRONT TO BACK)  
**30" / 760 MM**

LENGTH (FOLDED)  
**40" / 912 MM**

WIDTH (FOLDED)  
**20" / 507 MM**

DEPTH (FOLDED)  
**7.9" / 200 MM**

WEIGHT  
**37 LB / 14.7 KG**

LOAD CAPACITY  
**550 LB / 250 KG**

FERNO® RESERVES THE RIGHT TO CHANGE SPECIFICATIONS WITHOUT NOTICE.

## OPTIONAL ACCESSORIES

FOOTREST  
FLIP-UP HANDLES  
ARMRESTS  
HEADREST  
LINKAGE/CLOSE HANDLE

## STANDARD EQUIPMENT\*

WHEEL LOCKS  
TELESCOPING U-BAR  
RESTRAINTS