

FERNO EZ-GLIDE LBS Bariatric Stair Chair

Move patients of up to 318kg down stairs with ease –
and without risking operator injury



FOR MORE DETAILED
INFORMATION SEE OVER

Specifications

Height	Depth	Width	Weight	Load Limit
950 mm	850 mm	595 mm	31 kg	318 kg/50st/700lbs

FERNO – is a world leading innovator and manufacturer of patient transportation solutions.

Our product range is born of a closer understanding of our customers' needs and the practical demands of their work. Ferno's industry-leading **Bariatric EZ-Glide Evacuation Stair Chair** typifies our passion for top quality, highly functional design... patient transfer products that don't just do the job, they do it **BETTER**.

call +44 (0) 1274 851 999
email sales@ferno.co.uk
www ferno.co.uk

 **FERNO**[®]

FERNO EZ-GLIDE LBS Bariatric Stair Chair

The industry-leading patient handling of the EZ-Glide tracked chair...
now in a stronger, wider specification for larger patients

The EZ-GLIDE LBS combines all the industry-leading benefits of the tried and tested EZ-GLIDE tracked chair, but with an extended seat and strengthened fixed chair frame. This ensures that patients up to 318kg can be safely moved over a variety of surfaces, through confined areas and down stairways... all without risking injury to operators through strain and lifting.

This unique chair provides complete reassurance to the patient through its integrated track brake control system, which is so vital in potentially distressing evacuation scenarios.

- **Stronger and wider** – incorporating extended sides, reinforced footplate and a strengthened fixed frame, to seat patients with a greater physical mass
- **Uniquely easy-handling** – 25 mm track-to-ground clearance and smooth 100 mm front swivel wheels
- **Easily manoeuvrable in tight spaces** – 150 mm rear wheels for effortless handling even on most stair landings and other confined areas
- **Reassuring patient security** – with Biosafe 4 point harness and extending footplate complete with straps
- **Complete operator control** – unique adjustable braking system and rear carry handles provide further ease of handling
- **Feature rich** – including rear folding handles, extended footplate, IV pole and winch attachment, all included as standard
- **Improved comfort** – supplied with seat and back rest cushion, and headpad, headrest/shoulder support also available
- **Multi-purpose** – can be used as a in-vehicle patient seat when used with NMI Saferider® restraint system in conjunction with Pocket Lock® or Unwin Strap System®

Part Nos:

073132930 EZ-Glide LBS (Red) (Biosafe 4 point harness, with extended foot bar)

Replacement Parts and Accessories:

031402600 Replacement Seat cushion
FWE2110 Headrest/shoulder support
031402300 Dust cover
082197700 O2 cylinder holder

*For more information on NMI Pocket Lock and Saferider or Unwin Strap System, please contact our sales team



FERNO - Red EZ-Glide LBS FPC004 05.19

call +44 (0) 1274 851 999
email sales@ferno.co.uk
www ferno.co.uk



Ferno (UK) Limited, Stubbs Beck Lane, Cleckheaton, West Yorkshire, BD19 4TZ England.