

# FERNO **FRAC** Straps Pack & Immobilisers

The Emergency Services preferred, proven and efficient method of strapping or splinting damaged or fractured limbs



**FOR MORE DETAILED INFORMATION SEE OVER**

## Specifications

FRAC-Strap Pack	Length	Width	Depth	
1 x Wallet	470 mm	240 mm	110 mm	
1 x Padding Open	650 mm	100 mm	50 mm	
2 x Straps	1250 mm	130 mm		
2 x Straps	900 mm	75 mm		
1 x Strap	920 mm	130 mm		
FRAC-Immobilisers	Length	Height	Depth	Weight
5 Strap	750 mm	100 mm	150 mm	1.1 kg
4 Strap	600 mm	100 mm	150 mm	0.8 kg
2 Strap	350 mm	100 mm	150 mm	0.5 kg

**FERNO** – is a world leading innovator and manufacturer of patient transport solutions.

Our product range is born of a closer understanding of our customers' needs and the practical demands of their work. Ferno's industry-leading **FRAC STRAPS & IMMOBILISERS** typifies our passion for top quality, highly functional design... patient transfer products that don't just do the job, they do it **BETTER**.

**call** +44 (0) 1274 851 999  
**email** sales@ferno.co.uk  
**www** ferno.co.uk



# FERNO FRAC Straps Pack

**Frac-Straps** are easy to use and an efficient method of strapping breaks from collar bone to femur, helping prevent further injury before the casualty gets to A&E.

- **Extra comfort and tension support** - pack comes with 5 straps with elasticated area, 1 pad and storage wallet
- **Effective from application to treatment** - they are totally X-ray translucent
- **Easy to keep clean** - made from a tough, hygienic vinyl-coated nylon
- **Secure fitting** - straps have hardwearing VELCRO® fastenings



# FERNO FRAC Immobilisers

Effective box splints for suspected fractures at the scene of an accident. The extra padding also helps to prevent further damage to surrounding nerves and tissue during transportation.

- **Works on all shapes and size of patient** - three lengths accommodate all leg and arm size combinations
- **Effective from application to treatment** - they are totally X-ray translucent
- **Durable and reusable** - made from tough, hygienic vinyl-coated nylon for a long service life
- **Padded comfort and easy to apply** - open both ends U shape with large VELCRO® fastening straps help minimise trauma during application.



**Part Nos:**

031352501	Frac Straps Pack
031346901	Frac Immobiliser 2 Strap Child
031347201	Frac Immobiliser 4 Strap Short Leg or Arm
031347301	Frac Immobiliser 5 Strap Adult

FERNO - FRAC Straps & Immobilisers FPC005 10.15

**call** +44 (0) 1274 851 999  
**email** sales@ferno.co.uk  
**www** ferno.co.uk



Ferno (UK) Limited, Ferno House, Stubbs Beck Lane, Cleckheaton, West Yorkshire, BD19 4TZ England.